

Phone Adapter User Guide

Introduction

This user guide explains how the SoundGate is used for phone conversations.


The Phone adapter provides a wireless connection between your SoundGate and your home phone. It supports both incoming and outgoing phone calls. The SoundGate enables hands free operation.

First make sure the Phone adapter is correctly connected to both the phone land line and power. Please refer to the Phone Adapter Installation Guide.

The wireless phone application is possible within a range of approx. 30 feet (10 meters) from the Phone adapter.



SoundGate Overview



Phone calls can also be received when the keys are locked!

